



Oxivir Excel Foam

Cleaner and Disinfectant for non-invasive medical devices and nonporous surfaces

Description

Diversey Oxivir Excel Combined Cleaner & Disinfectant, powered by AHP Technology, provides you with the best alternative to your traditional disinfectants by delivering a full virucidal efficacy in just 30 seconds, while being very gentle to surfaces and with the lowest risk to your staff and customer. The active ingredient, hydrogen peroxide, breaks down into environmentally-friendly water and oxygen just minutes after use leaving no residues that might be harmful or change the appearance of the surface.

Key properties

- Full virucidal claim in just 30 sec
- Fast kill claims for a broad spectrum of viruses, bacteria and yeasticidal
- Superior Cleaning with high-performance surfactant system
- Fragrance free
- No rinsing required
- Safer for the environment, as well for the users, since the actives breaks down into water and oxygen as APEs, NPEs and VOCs

Benefits

- Broad efficacy provides the highest levels of disinfection. Extensively tested against European norms and effective against multiple pathogens of concern, carries a full virucidal claim including Polio, Adeno, and Murine norovirus, considered effective against all enveloped viruses such as: HIV, HBV, HCV, Influenza, Vaccinia and Coronaviruses.
- Safe for users no need PPE. No hazard classification
- Superior cleaning performance and 10x faster than your normal disinfectant
- AHP formulation is gentle on surfaces extending equipment life. No residues.
- Suitable for use in high risk areas in Retail, Education, Healthcare and other sectors
- Improve the disinfection compliance at 35%


Use instructions

Spray cleaning and disinfection:

1. Remove gross soils.
2. Apply to surfaces by either spraying the wiping cloth and wiping surfaces or applying directly to the surface by using a trigger sprayer.
3. After the contact time*, either wipe dry or allow to air dry.

Dry Wipes Disposable Wiping System:

**This dosage is according to optimal conditions, recommendations may vary, please consult with your Diversey representative for directions.*



© A.I.S.E.



Oxivir Excel Foam

1. Prepare Dry Wipes Disposable Wiping System according to instructions.
2. Remove gross soils.
3. Remove single wipe and wipe the surface.
4. After the contact time*, either wipe dry or allow to air dry.
5. Change the wipe when visibly soiled or moving to other cleaning area

**Ensure surface stays wet for duration of contact time.*

Technical data

Appearance: Clear, colourless liquid

Odour: Neutral

pH value (in use): 1.95

Relative density (20°C): 1.000

Flash point (°C): N/A

Shelf life: 720 days

The above data is typical of normal production and should not be taken as a specification.

Safe handling and storage information

Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet; sds.diverseym.com. Store in original closed containers away from extremes of temperature. Only for professional users / specialists.

Reacts with oxidisable substances. Shelf life of 2 years from date of manufacture.

Product compatibility

Important: Do not mix with other products. Do not use warm water. Do not use on water- or acid-sensitive surfaces. Under recommended conditions of use, Diversey Oxivir Excel Foam is suitable for use on most materials commonly encountered such as can be safely used and it is compatible with: stainless steel, chrome, vinyl, nylon fabrics, laminated surfaces, glass, rubber, hard and flexible plastics (polypropylene, polyurethane, polyethylene, PVC acrylic, fiberglass and polycarbonate)

Microbiological data

See EN table on last page

Environmental information

The surfactants used in this product are biodegradable in compliance with the requirement of EU Detergent Regulation, EC 648/2004.

Available pack sizes

Oxivir Excel Foam is available in 6x0.75L.



Oxivir Excel Foam

		Dirty conditions
	EN tests	Contact time
Bactericidal	EN 1276*	1 min
	EN 13727	30 sec
	EN 13697*	5 min
	EN 16615	1 min
Yeasticidal	EN1650*	5 min
	EN13624	5 min
	EN 13697*	5 min
	EN 16615	1 min
Virucidal	EN 14476 (Adenovirus, Poliovirus, Murine Norovirus)	30 sec
* Tests performed with dirty soiling for non-healthcare areas		