

## **Oxivir Plus**

### Oxivir Plus: Broad Spectrum Cleaner & Disinfectant

#### **Description**

**Diversey Oxivir Plus,** is a highly concentrated liquid cleaner disinfectant for the cleaning and disinfection of all water-resistant hard surfaces and some soft furnishings.

#### **Key properties**

- Hydrogen Peroxide based acidic formulation
- Effective against a wide spectrum of micro-organisms including bacteria, viruses, fungi and yeasts
- Effective cleaner
- Fragrance-free
- Based on patented technology; accelerated hydrogen peroxide

#### **Benefits**

- Cleans and disinfects in one step, reduction of working time
- Extensively tested against European norms and proven broad spectrum disinfection:
  Bactericidal, Virucidal, Fungicidal and Yeasticidal
- Excellent cleaning performance offering a streak free finish, descaling and lime limitation
- Suitable for use on water-resistant hard surfaces and some soft furnishings
- Good user safety profile
- · Good environmentally profile. Active breaks down into oxygen and water
- Suitable for use in all kind of application areas, eg in high risk areas in Healthcare

#### **Use instructions**

**Dosage:** User dilution 2-5% depending upon soil load and application need. Refer Microbiological data for efficacy of the product against multiple pathogens.

#### Application:

#### Spray cleaning and disinfection:

- 1. Remove loose dirt
- 2. Spray onto the surfaces to clean with clean cloth or pre-wetted with product.
- 3. If the surface is dry during cleaning, spray again to make the surface wet completely.
- 4. Allow to air dry.

#### Bucket cleaning and disinfetion:

- 1. Remove loose dirt
- 2. Apply solution with mop/cloth to clean and wet the surface completely.
- 3. Allow to air dry.

#### Important:

Do not use warm water. Do not mix with other products. Do not use on water or acid sensitive surfaces.







© A.I.S.E.



# Oxivir Plus

#### Technical data

Appearance : Clear yellow liquid

Relative density [20°C] : 1.029

pH-value in use  $: 1.8 \pm 0.5$  at 3.5% dilution

pH-value neat : 0.5

The above data is typical of normal production and should not be taken as a specification.

#### Safe handling and storage information

Full guidance on the handling and disposal of this product is provided in a separate Material Safety Data Sheet. Store in original container. Avoid extremes of temperature. Reacts with oxydisable substances.

#### **Product compatibility**

Under recommended conditions of use, Diversey Oxivir Plus is suitable for use on most materials commonly encountered.

#### **Environmental information**

The surfactants used in this product are ultimately biodegradable in compliance with the requirement of the EU Detergent Regulation, EC 648/2004.

#### Microbiological data

EN14476, Virucidal, Adenovirus, 2.0% 2 min, clean condition

EN14476, Virucidal, Poliovirus, 2.0% 2 min, clean condition

EN 14476, Virucidal, Human Adenovirus Murine Norovirüs, 2.5% 1 min, clean condition

EN14476, Virucidal, Human rotavirus, 1,0% 5 min, clean condition

EN14476, Virucidal, Vacciniavirus, 1.5% 5 min, dirty condition

EN 14476, Virucidal, Adenovirus, 2,5% 5min, dirty condition

EN14476, Virucidal, Poliovirus Adenovirus, 3,0% 5 min, dirty condition

EN1276, Bactericidal, Staphylococcus aureus Enterococcus hirae Escherichia coli Pseudomonas aeruginosa, 3,5% 30 sec; 3,0% 60sec, dirty condition

EN1276, Bactericidal, Staphylococcus aureus Enterococcus hirae Escherichia coli Pseudomonas aeruginosa, 2,0% 5min, dirty condition

EN13727, Bactericidal, Staphylococcus aureus Enterococcus hirae Pseudomonas aeruginosa, 1.0% 5 min, clean condition

EN13727, Bactericidal, Staphylococcus aureus Enterococcus hirae Pseudomonas aeruginosa, 2.5% 5 min, dirty condition

EN16615, Bactericidal, Staphylococcus aureus Enterococcus hirae Pseudomonas aeruginosa, 3.00% 5 min, dirty condition

EN 13727, Bactericidal, Staphylococcus aureus Pseudomonas aeruginosa Enterococcus hirae Acinetobacter baumannii, 0.25% 8 min, clean condition

EN13697, Bactericidal Yeasticidal, Staphylococcus aureus Enterococcus hirae Escherichia coli Pseudomonas aeruginosa Candida Albicans, Bactericidal: 1.5% 5 min Yeasticidal: 2.5% 5 min, clean condition

EN 13727, EN 13624 Bactericidal Yeasticidal Fungicidal, Escherichia coli: 1% 30 sec, Enterococcus hirae: 1% 30 sec,

Pseudomonas aeruginosa: 1% 30 sec, Staphylococcus aureus: 1% 30 sec, Candida albicans: 1% 5 min, Aspergillus niger: 1% 5 min

EN13624, Yeasticidal, Candida albicans, 2.0% 5 min, clean condition

EN13624, Yeasticidal, Candida Albicans, 3.0% 5 min, dirty condition

EN13697, Yeasticidal, Candida Albicans, 4.0% 5 min, dirty condition

EN14348, Mycobacteria, Mycobacterium terrae Mycobacterium avium, 2.5% 8 min, dirty condition

EN13727, Sporicidal, Clostridium difficile, 15% 10 min, clean condition

EN16615, Fungicidal, Aspergillus brasiliensis, 3.00% 15 min, dirty condition

EN13697, Fungicidal, Aspergillus brasiliensis, 3.5% 15 min, dirty condition