



Oxivir® Tb Wipes

Cleaning and disinfectant wipe for non-critical medical devices and surfaces (low level instrument grade disinfectant)

Oxivir Tb Wipes are intended for the cleaning and disinfecting of non-critical medical devices.



Based on proprietary AHP® technology to deliver fast, effective cleaning and disinfection.

Fast-acting

- Kills bacteria and viruses in just 1 minute
- Cleans and disinfects non-critical medical devices and healthcare environmental surfaces in 1 step

Effective

- Virucidal, Bactericidal, **Fungicidal and Tuberculocidal
- Oxivir Tb Wipes are proven to kill gram-positive and gram-negative bacteria and viruses including: MRSA, VRE, Enterococcus faecium, Listeria, Acinetobacter baumannii, Pseudomonas aeruginosa, CRE, Klebsiella pneumonia (CRKP), ESBL, MDR P. aeruginosa, PRSP, Adenovirus, Human Coronavirus 229E, Hepatitis B, HIV-1 (AIDS-virus), Influenza A (H1N1, H3N2), Norovirus, Poliovirus, Rhinovirus, Rotavirus, RSV, and Vaccinia Virus in 1 minute.
- Oxivir Tb Wipes kills Candida auris** in 5 minutes.
- **Kills SARS-CoV-2 (COVID-19 virus) in 15secs.**

Easy to use

- Ready-to-use pre-wetted wipes

Environmentally Responsible

- Active ingredient hydrogen peroxide, biodegrades into water and oxygen

Designed for

- Hospitals
- Aged Care facilities
- Child Care and Educational facilities
- Hospitality



AUST R: 338592



Oxivir® Tb Wipes

Cleaning and disinfectant wipe for non-critical medical devices and surfaces
(low level instrument grade disinfectant)

Efficacy Data:

Oxivir Tb Wipes are proven to kill Gram-positive and Gram-negative bacteria, enveloped and non- enveloped viruses, fungi and Tb.

Organism		Contact Time
Bacteria	Gram-positive	<ul style="list-style-type: none"> • Enterococcus faecium • Listeria • Methicillin-resistant Staphylococcus aureus (MRSA) • Vancomycin-resistant enterococci (VRE) • CA-MRSA (Genotype USA 300& 400) • Staphylococcus aureus
	Gram-negative	<ul style="list-style-type: none"> • Acinetobacter baumannii • Carbapenem-resistant Enterobacteriaceae (CRE) • Klebsiella pneumonia (CRKP) • MDR P. aeruginosa • MDR K. aerogenes • Enterobacter aerogenes • Escherichia coli • Escherichia coli 0157:H7 • ESBL • Klebsiella pneumoniae • Legionella pneumophila • Proteus vulgaris • PRSP • Pseudomonas aeruginosa • Salmonella enterica (choleraesuis) • Shigella Dysenteriae



Oxivir® Tb Wipes

Cleaning and disinfectant wipe for non-critical medical devices and surfaces
(low level instrument grade disinfectant)

Efficacy Data:

Oxivir Tb Wipes are proven to kill Gram-positive and Gram-negative bacteria, enveloped and non- enveloped viruses, fungi and Tb.

Organism		Contact Time
Viruses	Enveloped	<ul style="list-style-type: none"> • COVID-19 (SARS-CoV-2) • Human coronavirus 229E • MERS-CoV • Hepatitis B Virus (HBV) • Herpes simplex virus type1 & 2 • HIV-1 (AIDS virus) • Influenza A (H1N1 Swine flu, H3N2 Hong Kong) • Respiratory Syncytial Virus (RSV) • Vaccinia virus • Hepatitis C Virus (HCV)
	Non-Enveloped	<ul style="list-style-type: none"> • Adenovirus • Norovirus • Polio virus type 1 • Rhinovirus • Rotavirus
Fungi	<ul style="list-style-type: none"> • Candida auris • Candida albicans • Trichophyton Interdigitale • Trichophyton mentagrophytes 	<ul style="list-style-type: none"> • 5 mins • 3 mins • 5 mins • 10 mins
Tb	<ul style="list-style-type: none"> • Mycobacterium bovis (Tb) 	<ul style="list-style-type: none"> • 5 mins



Oxivir® Tb Wipes

How to Use Oxivir® Tb Wipes



READY TO USE PRE-WETTED CLEANING AND DISINFECTING WIPE

Use instructions

All surfaces must remain visibly wet for 1 minute.
To kill COVID-19, surfaces must remain visibly wet for 15 seconds.
For Candida auris, surfaces must remain visibly wet for 5 minutes.

Do not mix with other detergents or other chemicals. Not recommended for use on brass, copper or marble. Refer to Product Sales Sheet (PSS) for specific recommendations. Always follow medical device manufacturer's recommendations. Not intended for use on any invasive surgical instruments.
Do not flush wipe in a toilet. Dispose of contents/container/wipe in accordance with local waste management guidelines.

Technical data

Appearance	Pre-wetted wipes	Shelf Life	3 years
Relative density (20°C)	1.01		
Odour	Characteristic odour		
Active Ingredient	Hydrogen Peroxide		

The above data is typical of normal production and should not be taken as a specification.

Product	Pack size	Article code	Country
OxivirTb Wipes	12 x 160 wipes (15.5 x 17.3cm)	101104557	AU
Oxivir Tb Wipes-Large - Tub	4 x 160 wipes (28 x 30cm) Tub	100959102	AU
Oxivir Tb Wipes-Large –Refill	4 x 160 wipes (28 x 30cm) Refill for Tub	100959104	AU

Accelerated Hydrogen Peroxide® and Design are trademarks of Diversey, Inc.. All rights reserved.

Safe handling

Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet.

Storage information

Store in original container. Avoid extremes of temperature and humidity. Do not mix with other detergents or other chemicals. Store below 25°C.

Safety Reminder

Please make sure your employees read and understand the product label and Safety Data Sheet before using this product. The label contains directions for use, and both the label and SDS contain hazard warnings, precautionary statements and first aid procedures. SDSs are available online at www.diversey.com or by calling Toll free: Australia, 1800 647 779 or New Zealand, 0800 803 615.