



# Oxivir<sup>®</sup> Tb



## Hospital Grade Disinfectant

Oxivir Tb is a one step disinfectant, cleaner and odour neutraliser for cleaning and disinfecting most surfaces found in health, educational and recreational environments.

Oxivir Tb is proven kill gram-positive and gram-negative bacteria and viruses including: MRSA, VRE, Enterococcus faecium, Listeria, Acinetobacter baumannii, Pseudomonas aeruginosa, CRE, Klebsiella pneumonia (CRKP), ESBL, MDR P. aeruginosa, PRSP, Adenovirus, Human Coronavirus 229E, Hepatitis B, HIV-1 (AIDS-virus), Influenza A (H1N1, H3N2), Norovirus, Poliovirus, Rhinovirus, Rotavirus, RSV, and Vaccinia Virus in 1 minute.

**Kills SARS-CoV-2 (COVID-19 virus) in 15 seconds.**

### Effective

- One-Step Hospital Grade Disinfectant Cleaner based on Accelerated Hydrogen Peroxide (AHP<sup>®</sup>) Technology to deliver fast, broad spectrum disinfection with enhanced cleaning power.
- Virucidal, Bactericidal, \*\*Fungicidal, Tuberculocidal
- Kills Candida auris\*\* in 5 minutes.
- Sanitises soft surfaces such as curtains, cushions and carpets

### Easy to use

- Ready to use, non-rinsing formula
- One step cleaning and disinfection of hard environmental surfaces

### Safe

- Surface protection: non-corrosive, non-bleaching formula. Disinfect hard surfaces without creating destructive micro-cracks that harbour pathogens
- High safety profile – not classified as hazardous according to GHS requirements
- Oxivir Tb does not contain Methylisothiazolinone (MI), a preservative associated with systemic toxicity, dermatitis, skin corrosion and eye damage

### Environmentally Responsible

- No VOC's, no added fragrance
- Active ingredient, Hydrogen peroxide biodegrades into water and oxygen

### Designed for

- Healthcare (Hospitals, Nursing Homes, Medical Establishments)
- Education and Child Care Centres
- Hospitality, Offices and Retail



AUST R 165058



# Oxivir<sup>®</sup> Tb



## Hospital Grade Disinfectant

### Efficacy Data:

When used as directed, Oxivir Tb is highly effective against a wide variety of pathogenic micro-organisms including bacteria, antibiotic resistant bacteria, viruses, fungi and Tb.

|          |               | Organism  | Contact Time   |
|----------|---------------|---|--|
| Bacteria | Gram-positive | <ul style="list-style-type: none"> <li>• Enterococcus faecium</li> <li>• Listeria</li> <li>• Methicillin-resistant Staphylococcus aureus (MRSA)</li> <li>• Vancomycin-resistant enterococci (VRE)</li> <li>• CA-MRSA (Genotype USA 300&amp; 400)</li> <li>• Staphylococcus aureus</li> </ul>  | <ul style="list-style-type: none"> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> </ul>   |
|          | Gram-negative | <ul style="list-style-type: none"> <li>• Acinetobacter baumannii</li> <li>• Carbapenem-resistant Enterobacteriaceae (CRE)</li> <li>• Klebsiella pneumonia (CRKP)</li> <li>• MDR P. aeruginosa</li> <li>• MDR K. aerogenes</li> <li>• Enterobacter aerogenes</li> <li>• Escherichia coli</li> <li>• Escherichia coli O157:H7</li> <li>• ESBL</li> <li>• Klebsiella pneumoniae</li> <li>• Legionella pneumophila</li> <li>• Proteus vulgaris</li> <li>• PRSP</li> <li>• Pseudomonas aeruginosa</li> <li>• Salmonella enterica (choleraesuis)</li> <li>• Shigella Dysenteriae</li> </ul> | <ul style="list-style-type: none"> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 10 mins</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 10 mins</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> </ul> |
| Viruses  | Enveloped     | <ul style="list-style-type: none"> <li>• COVID-19 (SARS-CoV-2)</li> <li>• Human coronavirus 229E</li> <li>• MERS-CoV</li> <li>• Hepatitis B Virus (HBV)</li> <li>• Herpes simplex virus type 1 &amp; 2</li> <li>• HIV-1 (AIDS virus)</li> <li>• Influenza A (H1N1 Swine flu, H3N2 Hong Kong)</li> <li>• Respiratory Syncytial Virus (RSV)</li> <li>• Vaccinia virus</li> <li>• Hepatitis C Virus (HCV)</li> </ul>   | <ul style="list-style-type: none"> <li>• 15 secs</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 1 min</li> <li>• 5 mins</li> </ul>  |



# Oxivir<sup>®</sup> Tb



## Hospital Grade Disinfectant

### Use instructions

#### Ready to Use (RTU)

1. Apply Oxivir Tb undiluted to surfaces by spray, cloth or disposable wipe.
2. All surfaces must remain visibly wet for 1 minute for viruses and bacteria. To kill COVID-19, surfaces must remain visibly wet for 15 seconds. For Candida auris, surfaces must remain visibly wet for 5 minutes.
3. Allow to air dry.

A potable water rinse is required when disinfecting food contact surfaces. Do not mix with other detergents or other chemicals. Oxivir Tb can not be used on glasses, dishes and utensils or therapeutic devices. Not recommended for use on brass, copper or marble.

Hard surface disinfectant only. Not to be used on skin. Not intended to be used on medical devices or other therapeutic goods. Store below 25°C.

Oxivir Tb cleans, disinfects, deodorises and controls mould and mildew on hard non porous surfaces. Use daily on vinyl, plastic, glazed ceramic and porcelain, chrome, stainless steel, laminates, walls, ceilings, tables, chairs, countertops, telephones, other fixtures, toilet exteriors, patient rooms, operating rooms and shower rooms

### Technical data

|                         |                            |            |                      |
|-------------------------|----------------------------|------------|----------------------|
| Appearance              | Clear liquid               | Odour      | Characteristic odour |
| Specific Gravity (20°C) | 1.01                       | pH         | 3.0                  |
| Odour                   | Characteristic odour       | Shelf life | 2 years (RTU Liquid) |
| Active Ingredient       | 0.5% w/w Hydrogen Peroxide |            |                      |

*The above data is typical of normal production and should not be taken as a specification.*

| Product   | Pack size | Article code | Country |
|-----------|-----------|--------------|---------|
| Oxivir Tb | 12x946ml  | 5242201      | AU & NZ |

### Safe handling and storage information

Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet. Store in original container. Avoid extremes of temperature and humidity. Do not mix with other detergents or other chemicals. Store below 25°C.

#### Safety Reminder

Please make sure your employees read and understand the product label and Safety Data Sheet before using this product. The label contains directions for use, and both the label and SDS contain hazard warnings, precautionary statements and first aid procedures. SDSs are available online at [www.diversey.com](http://www.diversey.com) or by calling Toll free: Australia, 1800 647 779 or New Zealand, 0800 803 615.

For additional information, please contact us at:  
Diversey Australia Pty Ltd Toll free: 1800 647 779, Diversey New Zealand Ltd Toll free: 0800 803 615 [www.diversey.com](http://www.diversey.com). Legal disclaimer: The associated products or services are manufactured by or otherwise emanate from Diversey Companies and their subsidiaries and not from Sealed Air.

© 2020 Diversey, Inc. All Rights Reserved.