



TASKI®

Room Care



R2L

Hard Surface cleaner

Description

TASKI® Room Care R2L is a ready to use hard surface cleaner suitable for the cleaning of hard water resistant surfaces in hospitals, nursing homes, lodging and other establishments.

Features

- TASKI® Room Care R2L is registered as a Class I Medical Device, AUST L 367452.
- Non-abrasive.
- Non hazardous and non flammable.
- Pleasant, fresh fragrance.

Benefits

- Leaves surfaces hygienically clean.
- Quickly and gently removes oily and water-soluble residues.
- Suitable for most hard surfaces in ready to use squirt bottle.
- Colour and alphanumeric codes to prevent application mistakes.

Application areas

- Lodging
- Office Buildings
- Education
- Healthcare
- Retail





TASKI®

Room Care



R2L

Hard Surface cleaner

Directions for use

TASKI® Room Care R2L is ready to use.

1. Squirt a small amount of TASKI® Room Care R2L onto a damp cloth or paper towel and clean the surface.
2. Wipe or rinse with fresh water.
3. Buff with a dry cloth or paper towel for extra shine.

Technical data

Appearance	Clear green liquid	Odour	Fresh
pH (ready-to-use)	11.2	Specific gravity (20°C)	1.00
Dilution	Ready-to-use		

The above data is typical of normal production and should not be taken as a specification.

Product	Pack size	SKU Code	Availability
TASKI® Room Care R2L	6 x 750ml	4589756	AU

Safe handling and storage information

Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet. Store in original container. Avoid extremes of temperature and humidity.

Product compatibility

Under recommended conditions of use, TASKI® Room Care R2L is suitable for use on most washable hard surfaces. Test material compatibility in a small, inconspicuous place before use.

Safety Reminder

Please make sure your employees read and understand the product label and Safety Data Sheet before using this product. The label contains directions for use, and both the label and SDS contain hazard warnings, precautionary statements and first aid procedures. SDSs are available online at www.diversey.com or by calling Toll free: Australia, 1800 647 779 or New Zealand, 0800 803 615.

For additional information, please contact us at:

Diversey Australia Pty Ltd Toll free: 1800 647 779, Diversey New Zealand Ltd Toll free: 0800 803 615 www.diversey.com.

Legal disclaimer: The associated products or services are manufactured by or otherwise emanate from Diversey Companies and their subsidiaries and not from Sealed Air.