



Room Care R2-plus

Hard surface cleaner and disinfectant - concentrate

Description

Concentrated cleaner for the hygienic cleaning of all water resistant hard surfaces. To be used only with Divermite/DQFM dosing and dilution system.

Key properties

- Formulated for the hygienic cleaning of telephones and other hard, water resistant surfaces
- Slightly alkaline
- Non-abrasive
- Light, delicate and fresh fragrance
- Unique, patented ONT air neutralising technology alters malodour molecules

Benefits

- Leaves surfaces hygienically clean
- Quickly and gently removes oily and water-soluble residues
- Leaves a pleasant scent behind
- Neutralises malodour molecules, enhancing the effect of the perfume
- Colour and alphanumeric codes to prevent application mistakes



Directions for use

Use Room Care R2-plus at concentrations between 1-5 % v/v depending on the type and degree of soiling and application requirement.

When using with DQFM dilution dispenser unit, fill the Room Care R2-plus keyed spray bottle to the fill level mark directly from the DQFM dispenser. Securely screw-on the spray trigger. Proper dosage saves costs and minimises environmental impact.

Application:

Spray a small quantity onto a damp cloth/sponge pad and clean the surface. Let product act for the required time according to the microbiological data. Rinse with fresh water or wipe with a clean, damp cloth. Buff with a dry cloth for extra shine. For more detailed application, please refer to Room Care R2-plus Mini Guide.



©A.I.S.E.



Room Care R2-plus

Technical data

Appearance:	Clear, green liquid
pH value (neat):	≈ 11
pH value (in use):	>10
Relative density (20°C):	≈ 1.03

The above data is typical of normal production and should not be taken as a specification.

Safe handling and storage information

Full guidance on the handling and disposal of this product is provided in a separate Safety Data Sheet; sds.diverseyst.com.
Store in original closed containers away from extremes of temperature. Only for professional users / specialists.

Product compatibility

Do not use on water-sensitive surfaces, e.g. wood. Test material compatibility in a small, inconspicuous place before use.
Do not use with other dilution systems or methods. Do not use undiluted.

Microbiological data

Room Care R2-plus passed EN 1276 at 10% dilution in water of standardised hardness and in dirty conditions (0.3% bovine albumin) at a contact time of 1 minute.

Room Care R2-plus passed EN 13697 at 10% dilution in water of standardised hardness and in dirty conditions (0.3% bovine albumin) for referenced strains *Staphylococcus aureus*, *Enterococcus hirae*, *Escherichia coli* and *Pseudomonas aeruginosa* at a contact time of 5 minutes. Yeasticidal activity in 15 minutes when diluted at 10% for the referenced strain *Candida albicans*.

Environmental information

The surfactants used in this product are biodegradable in compliance with the requirement of EU Detergent Regulation, EC 648/2004.

Available pack sizes

Article code	Pack size	Pack type
100863768	2 × 1.5L	Pouch
JA102199	2 × 1.5L	Pouch